

ENDO



Ergonomic RA650 Front Protection Apron

<https://endo.id>


Ergonomic RA650

Front Protection Apron

Front protection in its pure form: Model RA650 offers full-body protection over the entire front section covering the required 60% of the girth.

Protection/Lead Equivalent Standard Collection

NovaLite, Leadfree or Standard lead



Front:

0.25 mm Pb, 0.35 mm Pb or 0.50 mm Pb

Back generally: 0.25 mm Pb

- Suitable for a wide range of applications with only frontal radiation exposure
- Ideal fit, maximum comfort for the wearer
- Possibility to slip off the apron
- Optionally "without protection" in the back, for a reduced weight of the apron
- Standard equipment includes shoulder pads
- Type testing and certification in accordance with (EU) 2016/425, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016



Sizes*, recommended for:

- XS Chest girth \leq 84 cm
- S Chest girth approx. 85 - 95 cm
- M Chest girth approx. 96 - 106 cm
- L Chest girth approx. 107 - 117 cm
- XL Chest girth approx. 118 - 128 cm
- XXL Chest girth approx. 129 - 139 cm
- 3XL Chest girth approx. 140 - 150 cm

Waist/hip girth \leq 91 cm

Waist/hip girth approx. 92 - 102 cm

Waist/hip girth approx. 103 - 112 cm

Waist/hip girth approx. 113 - 123 cm

Waist/hip girth approx. 124 - 134 cm

Waist/hip girth approx. 135 - 145 cm

Waist/hip girth approx. 146 - 156 cm

**The size depends on the girth of the chest, waist, or hip, depending on which of the three measurements is the largest.*

Closure System

ENDO uses only top quality, low-wear hook-and-loop fasteners.



Crossbands

Comfortable, wide crossbands allow the apron to be slid off easily by opening the hook-and-loop fasteners in the side area.

Outlast Temperature Regulation



Fast, Simple and Convenient

The apron can be slid off by simply opening the hook-and-loop fasteners in the side area, even through the sterile surgical gown.

Shoulder Relief

Thanks to the extra-wide, comfortable crossbands, the load is optimally distributed over the entire shoulder area.

Outer Material

Colours
(see front page)



Optional:
ComforTex® HPMF,
ComforTex® HPMF Hybrid

Sizes

From XS to 3XL
(see front page)



Lengths / Body height:
90 cm / 150 - 160 cm
100 cm / 160 - 170 cm
110 cm / 170 - 180 cm
120 cm / 175 - 190 cm
130 cm / 185 - 195 cm
140 cm / 195 - 205 cm

Optional



On request:
Outlast
Breast pocket
Embroidery
Made-to-measure
protection

Accessories


Bolero
Thyroid/sternum



ENDO

 Jalan Raya Menganti 14 Kedurus, Surabaya 60223
Jawa Timur, Indonesia

 62-31-7673636

 62-31-7673737

 0800 . 177 . ENDO (3636)

 info@endo.id

 <https://endo.id>