

ENDO

ENDO Oxygen Concentrator ENDO LFY-I-5A-01



Specification

Instrument Classification : Class II, Type B

Flow Rate : 10L/Minute, Float Flowmeter, Two Outlets

Oxygen Concentration : $\geq 93 \pm 3\%$ all flow rates

Noise : ≤ 55 dB (A)

Voltage / Frequency : 220V/50Hz, 110V/50Hz, 230V/50Hz

Power Consumption : 580W

Working Time : 24 hours non stop

Warranty : 1 year (starts after delivery)

Outlet Pressure : 0.03-0.07 Mpa

Working Temperature : 5°C - 40°C

Machine Size : 410 x 410 x 620mm

Packing Size : 495 x 475 x 730mm

Net Weight : 27 Kg

Gross Weight : 31 Kg



Supply 2 Patients



Triple Filtrations



24 Hours Non stop



Molecular Sieves



High Oxygen
Concentration



Low Concentration
Alarm



High Pressure
Alarm



Low Pressure
Alarm



Power Failure
Alarm