

ENDO

ENDO Oxygen Concentrator ENDO LFY-I-5F-11



Specification

Instrument Classification : Class II, Type B

Flow Rate : 5L/Minute, Float Flowmeter

Oxygen Concentration : $\geq 93 \pm 3$ all flow rates

Noise : ≤ 50 dB (A)

Voltage / Frequency : 220V/50Hz, 110V/50Hz, 230V/50Hz

Power Consumption : 400W

Working Time : 24 hours non stop

Warranty : 1 year (starts after delivery)

Outlet Pressure : 0.03-0.07 Mpa

Working Temperature : 5°C - 40°C

Machine Size : 520 x 550 x 200mm

Packing Size : 585 x 635 x 275mm

Net Weight : 21 Kg

Gross Weight : 22.5 Kg

Accessories : Nasal Vannula, Outer Filter, Humidifier Bottle

Optional : Remote Controller, Inner Filter, Nebulizer Kit,
Atomizer, Oxygen Purity Indicator



Atomizing Nozzle



Triple Filtrations



24 Hours Non stop



Molecular Sieves



High Oxygen
Concentration



Low Concentration
Alarm



High Pressure
Alarm



Low Pressure
Alarm



Power Failure
Alarm