




ENDO Shaver System

Arthroscopic Shaver



Features :

- Full touch LCD screen
- Easy-to-use interface, switch work mode effortlessly
- Two handpiece ports
- Working efficiency

Blades : oscilation speed range 500-3000RPM,

forward or reberse speed range 500-5000RPM

Burrs : forward or reverse speed range 6000-12000RPM



ENDO Shaver System

Blades and Burrs



Features :

- Single-use
- Blades and burrs available
- Different sizes and shapes
- Compatible with handpiece from Dyonics® and Livantec®

Ordering Information

	Blade / Bur	Compare to Smith & Nephew Dyonics®	Blade / Bur	Compare to Livantec®
	Teardrop Blade 4.2	7205310	Teardrop Blade 4.2	C9242
	Full Radius Blade 3.0	72201509	Full Radius Blade 3.0	C9941
	Full Radius Blade 4.2	7205306	Full Radius Blade 4.2	9247A/M9247A
	Inner serrated blade tooth blade 4.2	7210499	Inner serrated blade tooth blade 4.2	9263A/M9263A
	Double serrated tooth blade 3.0	72201513	Double serrated tooth blade 3.0	C9951
	Double serrated tooth blade 4.2	7205345	Double serrated tooth blade 4.2	C9254/MC9254
	Double serrated tooth plus blade 4.2	7210976	Double serrated tooth plus blade 4.2	M9299A
	Oval bur 5.9	7205327*/72200725	Oval bur 5.9	H9102/MH9102
	Round bur 5.9	7205325*/72200723	Round bur 5.9	H9112/MH9112