

ENDO



Model RA614

Thyroid & Sternum Protection

<https://endo.id>

Model RA614

Thyroid & Sternum

With Different Closures

Hook-and-loop fastener (RA614KL)

Standard size fits neck sizes of 32 - 45 cm.

Magnetic closure (RA614MA)


Standard size fits neck sizes of 32 - 45 cm.

Exchangeable hook-and-loop fastener (RA614AK)

Standard size fits neck sizes of 37 - 50 cm.


Protection/Lead Equivalent

Novalite, Leadfree
or Standard lead




As standard:
0.50 mm Pb

Outer Material



Optional:
ComforTex® HPMF,
ComforTex® HPMF Hybrid

Colours



Neck Sizes

Uniform
Depending on the model



Optional



On request:
Embroidery
Made-to-measure protection

Accessories/Replacement

Hygienic cover
(RA614HF)



Spare hook-and-loop fastener
(RA614AKV, adjustment
via snap buttons)

Thyroid and Sternum Protection

The thyroid is a radiation-sensitive organ of the human body and therefore in particular need of protection. Wearing an appropriate shield should therefore be a matter of course.

Our thyroid and sternum protector variants are characterised by a first-class cut and the best fit, taking into account all aspects of radiation protection.



Snap button connection

It allows for a secure adjustment of the thyroid and sternum protection on MAVIG X-ray protective clothing.



Accessories: Hygienic cover

To meet the hygienic requirements in everyday life, we recommend using our washable hygienic cover in the colour Titan made of the special ComforTex® MF material.

Model Variants RA614

The thyroid and sternum protector with snap button connection to the vest/coat/apron is suitable for all model series.

As a standard, it is manufactured with a firmly sewn-on hook-and-loop fastener. The other closure types, magnetic closure and exchangeable hook-and-loop fastener, are available as options. All variants provide enough room for adjustment.



Closure type-via magnet


The thyroid and sternum protection is also available with a practical magnetic closure (variant RA614MA).



Closure type – exchangeable

As a third alternative, an exchangeable closure can be chosen (RA614AK). Please note the different neck sizes here.

ENDO

 Jalan Raya Menganti 14 Kedurus, Surabaya 60223
Jawa Timur, Indonesia

 62-31-7673636

 62-31-7673737

 0800 . 177 . ENDO (3636)

 info@endo.id

 <https://endo.id>