

BERSAMA KITA ATASI COVID 19



Mari kita wujudkan kepedulian untuk para Pahlawan Kemanusiaan yang berada di garda terdepan melalui penyediaan produk APD yang ekonomis dan berkualitas

Alat Pelindung Diri

1 Set terdiri dari:

- Protective Clothes
 Materials:
 PP (Polypropylene fiber 50gr)
 PE (Polyethylene film) composite
- Face Shield
- Goggles: anti fog





Class 6

PT. ENDO Indonesia

Jl. Raya Menganti No. 14, Kedurus - Surabaya 60223





Tel: +62-31-767 3636

Fax: +62-31-767 3737



Email: info@endo.id URL: https://endo.id



Disposable Protective Clothing

Features

- Anti splash of liquid, body fluid and blood.
- Bacteria / Virus Blocking
- Waterproof
- Anti-Static

Material

PPPE "PP (50g) + PE", High grade for protection

Non-Sterile

• For medical staff, can be used in ICU, wards, outpatient, and other working environment which is exposed to high risk of infection

Packaging

- 30 pcs in a carton box
- Total weight is 8.6 kg,
- Carton box dimension is 0.438 * 0.308 * 0.508 (m)



Tight Cap Port

Prevent dust from entering.



Zip Design

Sealable storm flap for protection from contaminants.



Knitted Cuffs

For increased wearer comfort

No.	Height (cm)	Length, from top to bottom (mm)	Chest Circumference (mm)
S	160	1650	1200
	165	1690	1250
	170	1730	1300
XL	175	1780	1350
XXL	180	1810	1400
XXXL	185	1880	1450