



ENDO
INDONESIA

PRODUCT CATALOG

MEDICALS IMPORTER, DISTRIBUTOR, AND MANUFACTURER
Electromedic | Implant | Disposable | Paper

COMPANY PROFILE

ENDO INDONESIA was founded in 2006 as a company engaged in importation and distribution of Medical Equipment. Today our products include Electromedical Equipment, Reagents, Reusable and Disposable Medical Devices, among many others and accessories. And now we had grown into the manufacturing field of Medical Recording Chart Paper.

ENDO is a company with adequate legal and approval by the company's legal entity RI with number C-05336 HT.01.01. TH.2006. Supported by experienced technicians and proper equipments, we also provide repairs and calibration for electromedical equipment. ENDO exists to serve your various needs of Medical Equipments and highly concerned about customer satisfaction to promote long term working relationship.



TUV Nord ISO 9001 : 2008
Certified Company

 **Jl. Raya Menganti No. 14 Surabaya 60223,
Jawa Timur - Indonesia**

 **031 767 3636**

 **031 767 3737**

 **info@endo.id**

 **https://endo.id**

VISION

Becoming a valuable distributor of reputable manufacturers,
Continuously increasing employee's well-being,
Improving health care of Indonesian society
by providing good products at affordable prices.

MISSION

Consistent development and maintains integrity in daily routines,
while continuously providing best-value products and services.

VALUE

-  **Improving Healthcare**
Providing high quality products at affordable prices
-  **Networking**
Numerous sub-distributors nation-wide
-  **Qualified Personnels**
Experienced staffs in their respective fields
-  **Fast Delivery**
Most products are ready in stock

*Good Quality Products at
Affordable Price*



■ ENDO TRAUMA IMPLANTS

- *Plate & Screw Non Locking*
- *Plate & Screw Locking*



■ ENDO INTERLOCKING NAILS



ENDO Suction Pump 7A-23B



ENDO Suction Pump 7A-23B is a mobile suction apparatus developed from similar products and equipped with new generation oil-free negative pressure pump. Electric suction apparatus is applicable for suction of purulence and 'viscosity fluid. It's not applicable for other use and using by non-medical staff.

ENDO Suction Pump 7A-23B type of electric suction apparatus adopts completely plastic panel design, which makes it more fancy and fashionable. The device of intercepting oil for circumfluence can reduce the oil pollution. Easy to carry, low noise, high negative pressure and large flux. It can be widely applied in the surgical operations which need to absorb phlegm, thick liquid etc.

Technical data:

Power Voltage	: AC220V±22V,50Hz±1Hz
Max negative pressure	: ≥0.09MPa
Noise	: ≤65 dB (A)
Power	: ≤180VA
Pumping rate	: ≥40L/min
Reservoir capacity	: 2500mL/pc, 2pieces, glass
Vacuum Range	: 20kPa
Working System	: Continuous Operation
Packing Dimension	: 35.0×30.5×79.5cm, N.W:19.5KG



Large Flow (40 L)



Oil pump free

ENDO Portable Phlegm Suction Pump 7E-C



ENDO Portable Phlegm suction pump 7E-C is based on similar products at home and abroad is a new generation of oil free lubrication suction device, which is suitable to use by the patient who has difficulty in phlegm removal due to illness, coma and operation, as well as for aspirating such liquid as pus and blood during the clinical practice. It is the commonly applied medical device in the emergency room, operation room and for nursing in sickroom and home health care. ENDO Portable Phlegm suction pump 7E-C is designed with streamline shape, easy to carry, with heat release system, circular negative pressure meter, half hidden type liquid storage bottle and detachable power line.

Technical data:

Power supply	: AC220V±22V,50Hz±1Hz
Input Power	: 90VA
Max negative pressure value	: ≥0.075MPa
Pumping rate	: ≥15L/min
Noise	: ≤65dB (A)
Reservoir volume	: 1000mL, 1 piece
Working System	: Short Running
Packing Dimension	: 37.0x18.0x28.0 (mm), N.W:4.4KG



Medical Grade



Stepless Pressurize



Low Noise



Concealed Bottle

ENDO-NV8 Neonatal Ventilator

Specially designed for neonates, NV8 provides the most professional and diversified noninvasive nasal ventilation modes with advanced apnea wakeup function and automatic leakage compensation function.

Accurate and safe fresh gas delivered by iFlow Intelligent Closed-loop Control System to protect fragile new life; User-friendly operating system with easy clinical practice to facilitate medical staff; Integrated ventilation solutions providing SNIPPV/NIPPV, NCPAP and HFNC modes for safe ventilator weaning.

- ✓ User friendly design
- ✓ Safe and reliable
- ✓ Quick parameter setting, greatly reducing workload
- ✓ The world's leading ventilation modes



ENDO Vacuum Tubes offers such as blood collection system, swab ,specimen container ,six major series, high-quality security products for clinical laboratory, it has passed ISO13485, ISO9001 and CE approval, in full compliance with international and industry standards. Innovative production technology, product quality and stability.

ENDO vacuum blood collection system have good compatibility, can adapt the transmission device of various biochemical, immune, blood analyzer, using high quality imported additives products make the blood test is faster, higher accuracy.

- ✓ Vacuum Tube Sterile K3EDTA
- ✓ Vacuum Tube Sterile Gel & Clot Activator
- ✓ Vacuum Tube Sterile Lithium Heparin
- ✓ Vacuum Tube Sterile Plain
- ✓ Vacuum Tube Sterile Clot Activator
- ✓ Vacuum Tube Sterile Sodium Citrate



ENDO

Medical Recording Chart Paper




ECG Paper · Treadmill Paper · CTG / NST Paper
Spirometer Paper · EEG Paper · Patient Monitor Paper
Defibrillator Paper · Audiometer Paper · Urine Analyzer Paper
Hematology Analyzer Paper · Auto Refractometer Paper

Advantages


- Clean die cut
- Good ink adhesiveness (no ink residue on printhead)
- Micro perforation for easy tearing
- Precise gridline spacing
- Dark and vivid print out
- Accurate sensor mark
- Premium grade thermal papers
- Accurate precision
- Consistent quality and supply reliability
- Packing in dust free and hygienic environment
- More than 90 paper types, compatible for use on most major machine brands



ECG PAPER



TREADMILL PAPER



CTG/NST PAPER



THERMAL PAPER



CLEAR ULTRASOUND GEL

Good Quality Ultrasound Gel / USG Gel with affordable price.

Based on the same "can't be copied" formula of the World Standard; Ultrasound Transmission Gel, without fragrance or color. For diagnostic and therapeutic medical ultrasound procedures, like obstetrics, gynecology, digestion urinary and nervous systems, etc.

- No fragrance
- No color / transparent
- Acoustically correct for the broad range of frequencies used
- No formaldehyde

- Hypoallergenic, bacteriostatic, non-irritating
- 250ml / bottle (8.5 oz)
- 40 bottles / Carton Box
- Packing Size: 38cm X 27cm X 23cm (L x W x H)



DISPOSABLE ECG ELECTRODE

ENDO Disposable ECG Electrodes, which are used for various ECG tests. It uses an Ag/AgCl sensing element and Hydro-gel for adhesion. These are components for sensitive monitoring. Hydro-gel's main advantage is one of reducing movement artifacts by increasing the bond between the skin and the sensing element. The ECG electrode of ENDO ECG fields. An Ag/AgCl sensing element has the a best sensitivities and the solid adhesive hydro-gel has very low impedance.