




MEDICO ELECTRODES



Disposable ECG Electrodes

We Care For Your Heart !

HOLTER AND STRESS TEST ELECTRODES

A. Holter Electrodes



MLGLT-01, Liquid | MSGLT-01G, Solid
Size: 57 x 38 x 1 mm
Foam



MLGLT-03, Liquid | MSGLT-03G, Solid
Size: 50 x 48 x 1 mm
Foam



MLGLT-15, Liquid | MSGLT-15G, Solid
Size: Ø 55 x 1 mm
Foam

B. Stress Electrodes



MLGLT-02, Liquid | MSGLT-04, Solid
Size: 50 x 38 x 1 mm | 50 x 35 x 1 mm
Foam



MSGLT-08G, Solid
Size: 45 x 48 x 1 mm
Foam




MSGLT-09, Solid
Size: 51 x 33 x 1 mm
Foam



MLGLT-14, Liquid | MSGLT-14G, Solid
Size: Ø 50 x 1 mm
Foam

- Standard Packaging : 5 Electrodes/Card, 50 Electrodes/Pouch, 10 Pouches/Box, 12 Boxes/Carton.
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "h x w x t" format



For Illustration Purpose only

**ALL
PRODUCTS**



Sensor
Ag+AgCl

RESTING & MONITORING ELECTRODES

A. Diagnostic Resting Electrodes



MLGST-64, Liquid, | MSGST-64, Solid
Size: 44 x 32 x 1 mm
Foam, Banana Adaptor



MSGST-64M, Liquid, | MSGST-64M, Solid
Size: 44 x 28 x 1 mm
Foam, Banana Adaptor



TAB ELECTRODES
(See Overleaf)

B. Monitoring Electrodes



MSGLT-05MG, Solid
Size: 40 x 36 x 1 mm
Foam



MSGST-06, Solid | MSGST-06M, Solid
Size: 44 x 32 x 1 mm | 44 x 28 x 1 mm
Foam



MSGLT-15DC, Solid
Size: Ø 55 mm
Non-Woven Cloth




MSGLT-17T, Solid
Size: 49 x 48 mm
Clear Tape

- Standard Packaging : 5 Electrodes/Card, 50 Electrodes/Pouch, 10 Pouches/Box, 12 Boxes/Carton.
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "h x w x t" format

SUBSTRATE OPTIONS



For Illustration Purpose only

CONNECTOR OPTIONS



GEL OPTIONS



PEDIATRIC & NEO-NATAL ELECTRODES

A. Pediatric



MLGST-07, Liquid Foam
Size: 35 x 27 x 1 mm



MSGST-37, Solid Foam
Size: 35 x 30 x 1 mm



MSGST-27, Solid Foam
Size: 32 x 25 x 1 mm

- Standard Packaging : 5 Electrodes/Card, 50 Electrodes/Pouch, 10 Pouches/Box, 12 Boxes/Carton.

B. Neo-Natal - Prewired


Happy Trode



MPGBW - 30F, Hydrogel
Size: Ø 30 mm
Foam



MPRYB - 78, Hydrogel
Size: 22 x 22 mm
Non-Woven Cloth




MPRYB - 79RTT, Hydrogel
Size: 32 x 22 mm
Tape
Radiolucent



Happy Trode Plus⁺





MSRYB - 28AT, Hydrogel
Size: 28 x 20 mm
Non-Woven Cloth
Radiolucent



MSGBW - 31AT, Hydrogel
Size: 19 x 16 mm
Non-Woven Cloth
Radiolucent

- Standard Packaging : 3 Electrodes/Card, 3 Electrodes/Pouch, 100 Pouches/Box, 12 Boxes/Carton.
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "h x w x t" format



TAB ELECTRODES

A. Wet Gel Tab Electrodes



TL-3728
Size: 37 x 28 mm
Foam
Radiolucent

- Standard Packaging : 10 Electrodes/Card, 100 Electrodes/Pouch, 20 Pouches/Box, 6 Boxes/Carton.

B. Solid Gel Tab Electrodes



3426MED
Size: 34 x 26 mm
Paper / PET
Radiolucent



3423MED
Size: 34 x 23 mm
Paper / PET
Radiolucent





3420MED | 3421MED | 3521MED
Size : 34 x 20 | 34 x 21 | 35 x 21 mm
Paper / PET
Radiolucent



3413MED
Size: 34 x 13 mm
Paper / PET
Radiolucent

- Standard Packaging : 10 Electrodes/Card, 100 Electrodes/Pouch, 10 Pouches/Box, 10 Boxes/Carton.
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "hxwxh" format

1 HYGIENE	Prevents Nosocomial Infection: Effectively eliminates the possible risks of cross contamination related to electrodes use	
2 GOOD ADHESION	An effective adhesion throughout the ECG procedure	
3 TRACE QUALITY	Skin impedance is effectively and immediately reduced by using Tab Electrodes, giving instantaneous and reliable signal pick up	
4 PATIENT COMFORT	Safe and comfortable for all patient types	
5 EASE OF USE	Easy to place and remove, optimising the ECG Procedure and improving clinical workflows	
6 COST EFFECTIVE	Design is optimized to ensure cost effectiveness, saving valuable resources of the hospitals	OLD METHOD  NEW TECHNOLOGY  Quality Efficiency Cost

OFF-SET ELECTRODES

A. Stress Electrodes



OTL-4227-S42 : Ø 42 mm
OTL-4227-S50 : Ø 50 mm
OTL-5235-S60 : Ø 60 mm
Foam

- Standard Packaging : 1 Electrode/Card, 50 Electrodes/Pouch, 20 Pouches/Box, 6 Boxes/Carton

B. Resting Monitoring Electrodes



OTL-5235-S : Stainless Steel Stud
OTL-5235 : Tab Connector
OTL-5235-S4 : Banana Adaptor
Size : 52 x 35 mm
Foam



OTL-4227-S : Stainless Steel Stud
OTL-4227 : Tab Connector
OTL-4227-S4 : Banana Adaptor
Size : 42 x 27 mm
Foam

- Standard Packaging : 5 Electrodes/Card, 60 Electrodes/Pouch, 20 Pouches/Box, 6 Boxes/Carton
- Customization Options Available
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "h x w x t" format

1 Unique Off-Set connector
reduces motion artifacts

2 Comfortable foam backing with
enhanced adhesion ensures optimal
fixation throughout the procedure

3 Highly conductive wet gel
ensures reliable traces

4 Also available in
Radiolucent Connector

**ALL
PRODUCTS**



Sensor
Ag+AgCl

SPECIALIZED ELECTRODES

A. Tab & Snap Electrodes



TSL-4532
Size: 45 x 32 mm
Foam
Radiolucent

1

Radiolucent Applications

2

Supports both Tab Clip & Snap
Pinch Connector

3

Resting / Monitoring



For Illustration Purpose only

- Standard Packaging : 5 Electrodes/Card, 60 Electrodes/Pouch, 20 Pouches/Box, 6 Boxes/Carton

B. Repositionable Electrodes



REP4032RT, Hydrogel
Size: 40 x 32 mm
Non-Woven Cloth
Radiolucent

01

Radiolucent Connector

Allow "Re-Positioning" of the electrode
(Immediately after the first
use and on the same Patient)

02

03

Resting / Monitoring

For Illustration Purpose only

- Standard Packaging : 3 Electrodes/Card, 3 Electrodes/Pouch, 200 Pouches/Box, 12 Boxes/Carton

All Product Sizes are mentioned in "h x w x t" format

- Customization Options Available
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.



MEDICO ELECTRODES INTERNATIONAL LIMITED

Works: 142A/11&12, Noida Special Economic Zone, Noida 201305, India
Sales@MedicoElectrodes.com | www.MedicoElectrodes.com