



Always connected wireless,  
with extremely low battery consumption.

# Smart One<sup>®</sup>

Peak Flow and FEV1  
now available on your Smart Phone

*Simple and lightweight to carry on the go*  
Smart One<sup>®</sup> has been designed to satisfy the growing demand for wireless devices capable to connect to the new generation of Android and iOS Smart Phones and Tablets



*Ideal for monitoring respiratory illnesses and in the self management of Asthma, COPD, Lung Transplant care, Cystic Fibrosis and for use in Clinical Trials*






MIR SMART ONE



*Disposable turbine flowmeter  
Unrivaled in breath measurement.  
High accuracy and hygiene are 100% guaranteed.*



*MIR Smart One App includes an exclusive Incentive for adults and children based on both exhaled Flow and Volume that is helpful in improving subject's compliance*



*Graphic trend of Peak Flow or FEV1  
Symptoms selection and easy scoring system with the option of multiple tapping to identify the intensity*

