




SEESHEEN

HD image

Flexible Video Rhinolaryngoscope



EN Series


HD image



Zuhai Seesheen Medical Technology Co., Ltd.

EN Series


HD Image



High definition image shows more details of mucosal tissue.



Portable:
Compact design for ENT
Everywhere



CMOS+LED
Chip-on-tip Technology



Multi-function Buttons
Button ①: zoom in/out
Button ②: photo/video
Button ③: freeze



Application:

EN-1039: Rhinolaryngological examination and diagnosis in E.N.T.

EN-1249: Rhinolaryngological examination, diagnosis, and therapy in E.N.T.

Model	Insertion Tube	Channel	Working Length	Bending (Up/Down)	Field of View	Depth of Field
EN-1039	3.4mm	/	300mm	130°/130°	90°	3-50mm
EN-1249	4.8mm	2.0mm	365mm	130°/130°	90°	3-50mm

© Zhuhai Seesheen Medical Technology Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.
Add.4F,BLDG B, 6th Plant, Pingxi 10th Road, Nanping Science & Technology Park, Zhuhai 519060 China.
Tel: +86-755-23761871 Fax: +86-756-8676349 E-mail: intl@seesheen.com Web: www.seesheen.com

SEESHEEN

C-EN(HD)-202205-V1.1