



- Main Frequency : 400, 500kHz
- Power Source : AC120V or AC230V, 50/60Hz
- Protection Class : Class 1, Type CF
- Dimension (WxDxHmm) : 388x410x160
- Weight : 12.5kg

## Feature

- **Color Graphic User Interface**

Wide TFT LCD mount, clean and sharper image quality, unlike the existing electricity and the innovative new U.I. configuration operation.

- **Best Safety with Feedback System**

REM(Return Electrode Monitoring) function monitors the size of contact area between a patient and the pad. If the size is inappropriate, it automatically blocks the high-frequency current to minimize the danger of burning incidents.

- **7inch touch Screen System & key**

Touch screen system through the implementation easier and faster with a parameter setting, external key management can be the same as a touch screen.

- **Diversity of Feathre Bipolar**

Bipolar Cut, Bipolar Blend, Bipolar Forced to add features to implement more up-grade the optimal system.

- **Memory Function**

99 storage are available.

### Option 1 Endo Cut (Intra Cut)

Endo-cut is a monopolar electrosurgical treatment method and a cut which is carried out intervals (fractionated cut) with software-controlled alternations between cutting and coagulation stages.

## Output Specification

### Basic

	Mode	Output Power
Cut	Pure	450W/200Ω
	Blend 1	250W/200Ω
	Blend 2	200W/200Ω
	Blend 3	150W/200Ω
Coagulator	Contact	120W/200Ω
	Spray	100W/300Ω
	Argon	100W/300Ω
Bipolar Coagulation	Standard	100W/100Ω
	Auto Start	100W/100Ω
	Forced	80W/100Ω
Bipolar Cut	Cut	120W/200Ω
	Blend	100W/200Ω

### Option 1

	Mode	Output Power
Endo Cut	Papillotomy	150W/100Ω
	Polypectomy	150W/100Ω
Coagulation	Soft	80W/100Ω