

| TECHNICAL FEATURES | | | | |
|----------------------------|--|----------|----------|----------|
| DESCRIPTION | Reusable praecordial electrode for ECG | | | |
| REF. | F9008SSC | F9015SSC | F9009SSC | F9016SSC |
| METALLIC PART DIAMETER | 15mm | 20mm | 24mm | 30mm |
| BALL DIAMETER | 25mm | | 33mm | |
| ELECTRODE MATERIAL COATING | SILVER/SILVERCHLORIDE | | | |
| BALL MATERIAL | PLASTIC (MEDICAL GRADE PVC) | | | |
| CONNECTION | UNIVERSAL screw (4 mm) and clip | | | |
| COMPLIANCE TO STANDARDS | Directive 93/42/EEC. Class I device | | | |

| PACKAGING | |
|----------------|---|
| SALE PACKAGING | 6 pcs in open-close pouch |
| STORAGE | Temperature: 0°C ÷ +50°C – Relative humidity: 20% ÷ 80% . |