

TECHNICAL FEATURES			
REF.	F9023SSC	F9024SSC	F90240SSC
DESCRIPTION	Reusable peripheral clamp electrodes for ECG		
ELECTRODE MATERIAL COATING	Argento Clorurato (Ag/AGCl) / SSC		
CONNECTION	UNIVERSAL screw (4mm) and clip		
COMPLIANCE TO STANDARDS	Directive 93/42/EEC. Class I device		

PACKAGING	
SALE PACKAGING	4 pcs in polyethylene pouch
STORAGE	Temperature: 0°C ÷ +50°C – Relative humidity: 20% ÷ 80%