

의료용구 제조허가 제642호
식약청(품목허가) 제00-89호



VISCO

FOCUS-9

Medical Video Scope



NEW

Medical Video Scope

Millenium type


Standard type

VISCO

Feature

- The easiest program in the World!
- Unlimited Saving the images.
- Compare the Examination images with data.
- One touch saving, searching and comparing.
- Overlay board for High quality capture image.
- Print out 1 to 4 images immediately.

- 세계제일의 편리성체제 (Easy program)
- 무한대의 화상저장 (Millenium type)
- 진료화면과 Data화면 즉시비교 가능 (Millenium type)
- 화면비교, 저장, 검색 One touch (Millenium type)
- 전용 화상 Overlay board (Standard type)
- 즉석에서 1/4화면 Printing가능 (Standard type, Printer 추가)



Medical Video scope (Model: FOCUS-9)

Millenium type

- **Computer System**
 - 17" LCD monitor
 - Processor
 - Image database program
- **Videoscope**
 - main & handpiece
 - 12v 100w hallogen
- **System Trolley**

Standard type


- **Color Monitor 14"(CRT)**
- **ViedoScope**
 - 12v 100w hallogen
 - Overlay board (1or4 Image)
- **System Trolley**

Accessory

- **Otoscope 0° with lens**
- **Demascope(10X, 50X)**
- **Irisscope**
- **Dentalprobe anus scope**
- **OBGYN probe**

Option

- **LCD minitor 17"**
- **14"monitor (CRT)**



Technical Specification

Power supply	90-130V AC or 200-240V AC, 50/60H2
Image sensor	1/3" color CCD
Hprizontal resolution	460 TV Lines
Signal processing	Digital Signal Processor
AGC/White balance	Auto
Min illumination	2Lux
ELC	ON
Output signal	Composite:1EA or Video
Lamp	12V/100W Hallogen Lamp
Light intensity	Adjustable(0-100W)
Overlay Board	Change 1 to 4 Images Immediately

Millenium Type



Standard Type



•Specification & desiges subject to change without prior notice