


XTRUDER™

THE DEVICE FOR EASY AND FAST APPLICATION EVEN OF HIGH VISCOSITY RESINS AND BONE VOID FILLERS



XTRUDER is the new Tecres device for application of acrylic resins and bone void fillers prepared with the SHAKIT mixer.


ADVANTAGES:

- easy to fill

Thanks to the scraper ring, all you have to do is insert the XTRUDER syringe into SHAKIT and the product will flow into the syringe until it is full (maximum capacity: 10 ml).

- easy to apply

The specially designed shape and size of the ergonomic gigrigrip allow very easy extrusion of low, medium and even high viscosity resins.



TECRES® 
Advancing High Technology

www.tecres.it